


DONOR–ACCEPTOR OLIGOPYRROLES

Marcin Stepień

Faculty of Chemistry, University of Wrocław, Wrocław, Poland

Systematic tuning of optical bandgaps in organic molecules can be achieved by homologation (oligomerization) of linear π -conjugated motifs, by ring expansion of π -conjugated macrocycles, and by extension of fused ring systems in two dimensions. A complementary approach relies on combining donor and acceptor (D–A) moieties, with diverse recent applications in small-molecule and polymer chemistry. The D–A paradigm is particularly suitable for the development of tunable building blocks, which can be constructed by judicious merging of existing electron-deficient and electron-rich motifs. A simple and potentially productive design of such a hybrid structure is achieved by combining naphthalenemonoimide (NMI, red) and pyrrole (blue), as shown below [1,2]. We will discuss the application of such pyrroles as building blocks for the synthesis of diverse polycyclic aromatics, including macrocycles [2], nanographenoids [3,4], and bipyrroles [5,6], some of which are characterized by persistent helicene-like chirality. These systems reveal rich redox chemistry, spanning multiple oxidation levels, as well as tunable optical signatures extending into the near infrared.



- [1] H. Zhylitskaya, M. Stepień, *Org. Chem. Front.* **2018**, *5*, 2395–2414.
- [2] H. Zhylitskaya, J. Cybińska, P. Chmielewski, T. Lis, M. Stepień, *J. Am. Chem. Soc.* **2016**, *138*, 11390–11398.
- [3] M. Żyła-Karwowska, H. Zhylitskaya, J. Cybińska, T. Lis, P. J. Chmielewski, M. Stepień, *Angew. Chem. Int. Ed.* **2016**, *55*, 14658–14662.
- [4] M. Navakouski, H. Zhylitskaya, P. J. Chmielewski, T. Lis, J. Cybińska, M. Stepień, *Angew. Chem. Int. Ed.* **2019**, *58*, 4929–4933.
- [5] M. Żyła-Karwowska, L. Moshniaha, Y. Hong, H. Zhylitskaya, J. Cybińska, P. J. Chmielewski, T. Lis, D. Kim, M. Stepień, *Chem. – Eur. J.* **2018**, *24*, 7525–7530.
- [6] M. Żyła-Karwowska, L. Moshniaha, H. Zhylitskaya, M. Stepień, *J. Org. Chem.* **2018**, *83*, 5199–5209.